

JUMPING SPIDER

— Powered by Geyi Medical —

Single-use Endoscopy Solution



*See there,
Be there.*



NMPA





CONTENTS

- Single-use Ureteroscope
- Single-use Cystoscope
- Single-use Hysteroscope
- Single-use Bronchoscope
- Single-use Choledochoscope
- Portable Image Processor

Single-use Ureteroscope

As Slim as
6.3Fr




NMPA



PRODUCT FEATURE

Slim 6.3 Fr Distal

To support smooth ureteral access and efficient navigation.




275° Responsive Bidirectional Deflection

Designed to provide the flexibility needed to access challenging anatomy with greater control.



Y adapter included

packaged with each scope.



Handle Lever

- 1:1 coaxial steering, more smooth operation.
- Two shortcut keys, support user customization, can be preset with up to 15 functions.

240g Ergonomic Lightweight Handle

Designed for comfortable grip and intuitive scope manipulation.

3.6 Fr Working Channel

Supports compatible accessories while helping maintain irrigation performance.

Clear Digital Visualization

Digital imaging supports reliable observation throughout the procedure.



PRODUCT PARAMETER

Model	Outer Diameter	Working Channel	Working Length (mm)	Bending Angle	Depth of Field (mm)	Viewing Angle	Resolution (pixel)
GY-UR6.3	6.3 Fr	3.0 Fr	670	275° (up & down)	3-50	120°	160,000
GY-UR7.5	7.5 Fr	3.4 Fr					
GY-UR8.4	8.4 Fr	3.6 Fr					
GY-UR9.3	9.3 Fr	3.6 Fr					

Single-use Cystoscope

One Scope
Two Roles



EU-MDR




NMPA



PRODUCT FEATURE




220° Bidirectional Deflection
U220° / D220° deflection supports controlled navigation and systematic inspection of the bladder.



As Slim as 8.4 Fr
Designed for smooth transurethral access with a slim insertion profile.

120° / 140° Field of View
Wide viewing angle supports efficient inspection of the bladder.



Suction Connector

Rotation Part
90° L & R rotatable

One Scope, Two Roles
8.4 Fr, diagnostic purpose;
15.6 Fr, therapeutic purpose.

Designed for Comfortable Handling
Ergonomic handle design supports intuitive and responsive scope manipulation.

Large Working Channel
Up to 6 Fr, provides room for compatible urological instruments during interventional procedures.

PRODUCT PARAMETER


Model	Outer Diameter	Working Channel	Working Length (mm)	Bending Angle	Depth of Field (mm)	Viewing Angle	Resolution (pixel)
GY-CYM084WJQB145 aka:GY-CYM084	8.4 Fr	3.6 Fr	400	220° (up & down)	3-50	120°	160,000
GY-CYN156WJQA134 aka:GY-CYN156	15.6 Fr	6.0 Fr				140°	1,000,000

Single-use **Hysteroscope**

No Dilation

No Anesthesia

See & Treat



EU-MDR | **CE**

NMPA



PRODUCT FEATURE

4 Types of Bending Angles

0°/12°/22°/30°, it makes operation more convenient.



Dual Shortcut Keys

Support user customization, can be preset with up to 15 functions.

Type-C Connector

Suction Tube

Irrigation Tube

Working Channel

Compatible with 5 Fr and 6 Fr instruments.

Ergonomic Handle

With anti-slip surface treatment.

Dual-channel design

Flow regulator allows for free adjustment of the flow rate; squeezing the tube to its minimum diameter prevents back flow of the liquid, making it easy to operate.

HD Image

- Integrated distal CMOS with dual LED lights delivers clear images.
- Resolution: 160,000 pixels / 1,000,000 pixels, delivers excellent color fidelity.



PRODUCT PARAMETER

Model	Outer Diameter (mm)	Working Channel (mm)	Working Length (mm)	Bending Angle	Depth of Field (mm)	Viewing Angle	Resolution (pixel)
GYUT-A 12.6-0	4.2	2.2	220	0°	3-50	120°	160, 000
GYUT-A 12.6-12				12°			
GYUT-B 15-0	5.0	2.0		0°	3-100	140°	1, 000, 000
GYUT-B 15-12				12°			

Single-use **Bronchoscope**

Airway, Always Ready



NMPA



PRODUCT FEATURE

Handle Lever

Uniform damping on handle lever, instant stop at any position.

Type-C Connector

Suction Button

Suction Connector

Custom Button

Ergonomic Handle

With anti-slip surface treatment.

Working Channel

Direct-access instrument channel enables convenient operation.

Rotation Part

90° L & R rotatable

High Sensitivity CMOS+LED Lights


HD Image

Megapixel clear field of view, assists clinical discrimination and observation of bronchial airway tissues.



220°

Bidirectional Bending




Flexible bronchoscope with precise deflection; multiple specifications to meet diagnosis & treatment demands of different patient groups.

PRODUCT PARAMETER

Model	Outer Diameter	Working Channel	Working Length (mm)	Bending Angle	Depth of Field (mm)	Viewing Angle	Resolution (pixel)
GY-BRM075A1H	7.5 Fr	3.6 Fr	600	220° (up & down)	3-50	120°	160,000
GY-BRM084B1GL	8.4 Fr	3.6 Fr	600				
GY-BRM120B1GL	12.0 Fr	6.0 Fr	600				
GY-BRN147B1GL	14.7 Fr	6.6 Fr	600		3-100	140°	1,000,000
GY-BRN147B1GH	14.7 Fr	6.6 Fr	430				

Single-use **Choledochoscope**

Into the Duct, Into View



NMPA



PRODUCT FEATURE

Handle Lever

Uniform damping on handle lever, instant stop at any position.

Type-C Connector

Ergonomic Handle

With anti-slip surface treatment.

Rotation Part

90° L & R rotatable

Suction Button

Custom Button

Working Channel

Direct-access instrument channel enables convenient operation.

High Sensitivity CMOS+LED Lights


HD Image

Megapixel clear field of view, assists clinical discrimination and observation of bronchial airway tissues.



220°

Bidirectional Bending



Flexible choledochoscope with precise deflection; multiple specifications to meet diagnosis & treatment demands of different patient groups.

PRODUCT PARAMETER

Model	Outer Diameter	Working Channel	Working Length (mm)	Bending Angle	Depth of Field (mm)	Viewing Angle	Resolution (pixel)
GY-CTM075A1L	7.5 Fr	3.6 Fr	600	220° (up & down)	3-50	120°	160,000
GY-CTM084B1L	8.4 Fr	3.6 Fr	600				
GY-CTM120B1GL	12.0 Fr	6.0 Fr	600				
GY-CTM120B1GH	12.0 Fr	6.0 Fr	400		3-100	140°	1,000,000
GY-CTN147B1GL	14.7 Fr	6.6 Fr	600				
GY-CTN147B1GH	14.7 Fr	6.6 Fr	400				

Portable Image Processor

One Tablet, Multiple Endoscopy Solutions.



NMPA



PRODUCT FEATURE

Ureteroscope
Single-use

Cystoscope
Single-use

Choledochoscope
Single-use

Bronchoscope
Single-use

Hysteroscope
Single-use

12.1" Touch Screen
Intuitive touch operation, 1280 × 800 Display / 1920 × 1080 Extended Output.


Compatible With Various Endoscopes

Type-C Connector

Up to 8 Hours Battery
Portable and flexible operation.

32 GB Built-in Storage
Convenient image & video export

Special Features



PRODUCT PARAMETER

Model	Video Output	Display Resolution	Built-in Battery
GY-AIO-121	HDMI1.4	1280×800 (720P)	Excluded
GY-AIO-121(A)	1920×1080 (1080P)		Included

GLOBAL DISTRIBUTORS WANTED



ZHEJIANG GEYI MEDICAL INSTRUMENT CO., LTD.

📍 No. 8, Wanqitaibao Road, Yangxi Subdistrict, Jiande City, Hangzhou, Zhejiang Province, China

✉ info@geyi-medical.com

🌐 www.geyi-medical.com

🌐 www.jumpingspider-med.com