

JUMPING SPIDER


— Powered by Geyi Medical —

Single-use Hysteroscope

No Dilation

No Anesthesia

See & Treat



EU-MDR | C E

NMPA



PRODUCT FEATURE

4 Types of Bending Angles

0°/12°/22°/30°, it makes operation more convenient.



Dual Shortcut Keys

Support user customization, can be preset with up to 15 functions.

Type-C Connector

Suction Tube

Irrigation Tube

Working Channel

Compatible with 5 Fr and 6 Fr instruments.

Ergonomic Handle

With anti-slip surface treatment.

Dual-channel design

Flow regulator allows for free adjustment of the flow rate; squeezing the tube to its minimum diameter prevents back flow of the liquid, making it easy to operate.

HD Image

- Integrated distal CMOS with dual LED lights delivers clear images.
- Resolution : 160,000 pixels / 1,000,000 pixels, delivers excellent color fidelity.



PRODUCT PARAMETER

Model	Outer Diameter (mm)	Working Channel (mm)	Working Length (mm)	Bending Angle	Depth of Field (mm)	Viewing Angle	Resolution (pixel)
GYUT-A 12.6-0	4.2	2.2	220	0°	3-50	120°	160,000
GYUT-A 12.6-12				12°			
GYUT-B 15-0	5.0	2.0		0°	3-100	140°	1,000,000
GYUT-B 15-12							

ZHEJIANG GEYI MEDICAL INSTRUMENT CO., LTD.

✉ info@geyi-medical.com

🌐 www.geyi-medical.com

🌐 www.jumpingspider-med.com